
Taema

Manuel d'utilisation	(F)
Handleiding voor	(NL)
Gebrauchsanweisung	(D)
User's manual	(GB)
Manual del usuario	(E)
Manual de utilização	(P)
Manuale di utilizzazione	(I)

FREELOX

SOMMAIRE

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CONTENTS

GENERAL SAFETY INSTRUCTIONS	1	5. In case of problem	7
STATIONARY & PORTABLE DIAGRAMs	2	6. Symbols - Abbreviations	8
I. INTRODUCTION	3	7. Standards	8
II. DESCRIPTION AND USE	3	8. Technical safeties	8
1. Description of stationary canister (Figure 1)	3	9. Material in direct or indirect contact with the patient	8
2. Description of portable canister (Figures 2a & 2b)	3	10. Operating principle	8
3. Description of the level indicator (Figure 3)	3	IV. APPENDIX	9
4. Use of stationary and portable canisters	4	Installation checklist / After maintenance checklist (before each patient installation)	9
5. Filling the portable canister	4	Starting checklist (before each use)	9
6. Battery replacement	5	Maintenance sheet	10
7. Disposal of waste	5		
8. Disposal of the apparatus	5		
III. WHAT YOU NEED TO KNOW	6		
1. Recommended accessories	6		
2. Cleaning and disinfection	6		
3. Travel	6		
4. What to do in case of an accident	7		

GENERAL SAFETY INSTRUCTIONS

USING OXYGEN

The following guidelines must absolutely be adhered to when in the presence of oxygen:

- Do not leave the stationary canister connected to the portable canister,
- Close the flow rate selector knob whenever the unit is not in use,
- Place the system in a ventilated area,
- Canisters must absolutely be kept in a vertical position.

Oxygen is not a flammable gas, but it accelerates the combustion of materials. To avoid any risk of fire, the **FREELOX** should be placed more than 1.5m from any smoker(s), flames, electrical equipment, combustible products such as: oils, greases, solvents, aerosol vaporizers, etc.

SAFETY WITH COLD SUBSTANCES

Liquid oxygen is extremely cold (-183°C). The parts of the tanks that have been in contact with liquid oxygen (in particular during filling of tanks) can result in a burn in case of contact with the skin.

To avoid burns, never touch cold or frosted parts, avoid any spilling and pouring of liquid by absolutely keeping tanks in the vertical position.


Use the **FREELOX** system only for oxygen therapy. Oxygen therapy should only be done according to a medical prescription, adhering to the prescribed flow rates.

In case of long period of non use or storage, remove the batteries.

Only persons having completely read and understood this manual are authorized to manipulate and use the FREELOX system.

Check the maximal distance and the buzzer volume in function of environmental noise in order to keep alarm audible.

This unit is in accordance with the directive EEC 93/42 (medical devices), but can be affected by the nearby of the following equipment :
cellular phones, Citizen Band, micro-wave ovens, and more generally by electromagnetic interferences exceeding the levels specified in CEI 60 601-1-2.

I. INTRODUCTION

The **FREELOX** autonomous oxygen therapy system delivers medical oxygen at a flow rate which is adjusted according to a prescription. It consists of an assembly made up of two parts:

A stationary canister whose capacity is 20, 32 or 44 liters which is mounted on a base with wheels to allow mobility. This stationary canister makes it possible both to supply the flow of oxygen of the prescription and fill a portable canister.

A portable canister whose capacity is 0.5 or 1.2 liters which gives you total mobility with a substantial oxygen autonomy (for example, seven hours of autonomy at a flow rate of 2 liters of oxygen per minute with the 1.2 liter portable canister).

In case of dependence of the patient with oxygen (24H/24 treatment), an emergency supplying procedure must be provided by the distributor.

- Should the **FREELOX** system appear not to function correctly under normal conditions of use, the distributor (whose address must be given on the last page of the current manual) is alone authorized to troubleshoot and repair the system. Any attempt at repair by the user can be harmful to his/her health or safety.

II. DESCRIPTION AND USE

1. Description of stationary canister (Figure 1)

- | | |
|---|---|
| ① | Level indicator |
| ② | Battery compartment |
| ③ | Filler opening |
| ④ | Flow rate selector knob |
| ⑤ | Humidifier oxygen outlet |
| ⑥ | Unlocking knob |
| ⑦ | Humidifier |
| ⑧ | Water collecting bottle |
| ⑨ | 20, 32 or 44 liters stationary canister |
| ⑩ | Base with wheels |

2. Description of portable canister (Figures 2a & 2b)

- | | |
|---|------------------------------------|
| a | Level indicator |
| b | Oxygen outlet |
| c | Flow rate selector knob |
| d | Battery compartment |
| e | Vent lever |
| f | Water collection bottle (or cover) |

3. Description of the level indicator (Stationary and portable canister) (Figure 3)

- | | |
|-----|--|
| I | Display |
| II | Sectors (Liquid O ₂ level indication) |
| III | AIR LIQUIDE I.D. |
| IV | Audible alarm cutoff pushbutton |
| V | Charge indicator (battery) |

4. Use of stationary and portable canisters


1) Check level indicator (1 & a) to see if the canister contains enough oxygen. If not, ask the distributor to come and refill the stationary canister.

Comments:

- A low level is programmed by the distributor.
- When the level drops below the minimum threshold, display unit blinks and a BEEP every minute is audible.
- To stop the audible alarm, press the key .
- After 5 minutes, the audible alarm is automatically cut off in order to protect the battery lifetime.


- An indicator shows the state of charge of the level indicator battery : .

- 1 : continuous display : charge OK
 2 : blinking display: replace battery
 3 : no indicator display: setting mode
 4 : no display: battery is out of order.


2) Check the water collecting bottle (8 & f) is clean and empty the put it into place.

3) Connect fitting of humidifier (12 : stationary canister only) to oxygen outlet (11).


4) Unscrew bottle (13) from humidifier and fill it with water to the fill line. (Cf the humidifier instructions).

5) Screw bottle under its lid. Squeeze enough the cover to block any leak.

6) Connect your oxygen tubing (15) (or connect directly the nasal cannula) to humidifier outlet nipple (5). The length of the tube between the **FREELOX** stationary canister must be inferior to **15 meters** to provide a sufficient oxygen flow rate. For the portable canister, the cannula must be directly connected to the outlet nozzle (b).


directly connected to the outlet nozzle (b).

7) Turn flow rate selector knob (4 & c) clockwise to the prescribed flow rate (in liters/minute up to 7 l/mn).


Comments:

- the selector contains a stop which limits flow to the prescribed rate. Do not force the knob beyond this limit.
- be sure that the selector knob is correctly positioned. The **FREELOX** delivers no oxygen flow when the selector knob is between two values.

8) Check that there is a steady flow of bubbles in humidifier (when used).

9) Check that oxygen is flowing out of the inhalation device (nasal cannula or other device) by placing a finger in front of the outlet opening(s) (or by placing the opening(s) on the surface of water. In latter case, surface movements materializes oxygen flow.


10) Adjust the inhalation device (nasal cannula or other device not shown).

11) At the end of use, turn off the unit by setting the flow rate selector knob to 0.

Note:

Empty water collecting bottle (8 : stationary canister) when it is half full. Clean it with soapy water, then dry it before use it again. The piece of felt in the condensation collection tray (portable canister) must be dried after each use.

5. Filling the portable canister

It is recommended to fill the portable canister just before use. Because oxygen in portable canister evaporates, even not in use.


1) Check that the level in stationary canister is greater than 1/4.

2) Remove all traces of humidity on filler opening (3) of the **FREELOX** assembly using a clean lint-free cotton cloth.


3) Set flow rate selector (c) of portable canister to 0.

Remove the piece of felt placed in the water collection cover (portable canister), then dry it after each use.

4) Remove water collection bottle (8 : stationary canister) or condensation collection tray (f : portable canister).


5) Align and lower the portable canister vertically onto opening (3) in the cover of the stationary canister.

6) Lock the transfer device by pressing vertically on the upper part of the portable canister in the axis of the fittings .

7) Lower vent lever (e) of the portable to 45°; at this moment, a dull hissing sound due to the escape of oxygen gas is heard.

8) The portable canister is full after approximately one minute, when the hissing due to the escape of gas changes and a dense white vapour escapes from the cover of the stationary canister.

9) Raise vent lever up to its original position.

10) Disengage the portable canister by simultaneously pulling upwards on the handle and pushing unlocking button (6).

11) The white vapour escape must stop in a while.

comments:

To avoid malfunction during filling of the portable, remove all humidity trace on the filling openings which may form ice blocking stationary filling outlet in open position.

In case fittings are stuck due to frost, **do not force** ; wait a few minutes for the frosted parts to warm up and disengage the portable canister as soon as the frost has melted.

In case of liquid oxygen leakage after disconnection of portable canister, immediately reconnect the portable canister to the stationary canister, then separate them again. If the leak persists, reconnect the canisters again and ask the distributor for assistance.


12) Check fill level (a) with electronic level indicator.

13) Reinstall condensation collection tray and the piece of felt inside.

14) At the end of use, turn off the unit by setting flow rate selector knob (c) to 0.


Note : After each use, leave the portable canister in a vertical position (even when empty).

6. Battery replacement

1) Swing open battery compartment lid (2 & d) using a flat screwdriver.


2) Lift out battery and its connector.

3) Disconnect battery and replace it with a new one.

4) Reinstall battery and connector into the battery compartment.

5) Reinstall battery compartment lid until it clicks into place.

6) Check battery charge state by verifying that battery symbol  is displayed steadily.

Note:

The stationary canister level indicator is powered by a type 6LF22 CEI 9 volt alkaline battery.

7. Disposal of waste

All waste arising from the use of the **FREELOX** (hoses, filters, batteries,) must be disposed of in accordance with the appropriate hospital procedures.

8. Disposal of the apparatus

With the aim of preserving the environment, all equipment disposal must be carried out in accordance with the appropriate hospital procedures.

In order to respect the traçability imposed by the **CE** marking, it is obligatory to communicate the serial number of the eliminated device to the **Taema** technical service.

III. WHAT YOU NEED TO KNOW

1. Recommended accessories

The accessories used with **FREELOX** must be :

- oxygen compatible,
 - bio-compatible,
 - conforming to the general requirements of the standard EN 60601-1 and the directive CEE 93/42, and must not be antistatic or conductor of electricity.
- Hoses, cannula, masks, ... must be specific for oxygen-therapy.

The maximal length of the hose between the patient and the stationary canister **FREELOX** must be inferior to **15 meters** to provide a sufficient oxygen flow rate.

The use of accessories that do not conform to the requirements above absolves the manufacturer from responsibility in the event of an incident.

Accessories given reference numbers by Taema or included in the batch of accessories supplied with the apparatus do conform to these requirements.

2. Cleaning and disinfection

2.1. Cleaning

Cleaning is limited to the outside parts of the **FREELOX** system.

Important

The use of alcohol, acetone, solvents or any other flammable product **is prohibited**.

- For the plastic parts, clean with a cloth or a sponge moistened, or with cloth impregnated with aqueous alcohol solution (do not use abrasive powder).

Using an aqueous alcohol solution imposes that canisters are not operating (selector knob on position 0).

- for the stainless steel stationary canister, use of current domestic products (ammoniated cream) used with a not abrasive sponge is suitable Rinse, then leave to dry before use.

2.2. Disinfection

- By the patient :

Disinfection concerns only the external accessories (humidifier, oxygen inhalation device).

Concerning inhalation device, see manufacturer's recommendations.

Steadily empty the humidifier bottle and the water collecting bottle, clean them with soapy water, rinse and leave them to dry.

Note : in case of long period of non use or storage, keep humidifier bottle and water collecting bottle **empty and dry**.

- For each new patient :

Humidifier must be sterilized if possible (case of **NEBAL** fitted with **FREELOX**) or replaced.

The piece of felt (portable) must be replaced by a new one. The water collecting bottles (plastic bottle for the stationary & cover for the portable) must be disinfected by immersion in a disinfecting solution.

3. Travel

3.1. Autonomy

The portable canister allows true autonomy for travel for several hours (see table below).

Autonomy is given for a 24/24h use in optimum filling conditions.

Autonomy (portable and main 100% full)					
flow rate l/min	Portable 1,2l	Portable 0,5l	Stationary canisters		
			44l	32l	20l
0,5	22h	10h	50d	36d	22d
1	10h	5h	25d	18d	11d
1,5	8h	3h	17d	12d	7d
2	6h	2h30	12d	8d	5d
3	4h	2h	8d	6d	3d
4	3h	-	6d	4d	2.5d
5	2h	-	5d	3d	2d

3.2. Recommendations

When not in use, the portable canister empties naturally by evaporation. It is recommended to use the portable canister after travelling with it in order not to waste oxygen.

The condensation collection tray is equipped with a piece of felt which must be dried after each use.

Do not carry the portable canister under clothes (jacket, coat, etc.).

3.3. Travelling by car

For car travel, even for very short distance, only the portable canister is authorized.

Place the portable canister in the vertical position, suspended from a headrest or strapped to a seatback with a belt.

Never put the portable canister in the trunk or in an enclosed area.

Keep a window open and do not allow anyone to smoke.

4. What to do in case of an accident

4.1. In case a canister is overturned

Stand the canister upright immediately.

If the canister has been overturned for several minutes, you must ventilate the premises, for example by opening a window.

Do not induce any flames or sparks,

Leave the premises at least 1/2 h

If the portable canister is tipped over in a motorized vehicle:

- switch off the engine,
- ask all passengers to leave the vehicle,,
- Thoroughly ventilate the vehicle before continuing travel

4.2. In case of liquid oxygen spray

In eyes:

- Wash eye(s) copiously with water for at least 15 minutes.
- Call a physician.

On skin :

- Do not rub. Remove or loosen clothing if necessary. Defrost affected parts by moderate warming (with warm water, if possible, otherwise by placing them against another warm part of the body).
- Call a physician.

5. In case of problem

Observations	Probable causes	Remedies
• Oxygen flow too weak or zero.	1 Flow rate selector knob badly positioned 2 Canister is empty 3 Defective tube(s) 4 Inhalation device defective (canula, ...) 5 Humidifier (or outlet nipple) obstructed or leaking	1 Turn the knob to the correct position. 2 Check level. If canister is empty : - for portable canister, fill it, - for stationary canister, call your distributor. 3 Disconnect the tube. If flow is restored, check the tube (leak, kink or obstruction). Replace it if necessary. 4 Disconnect inhalation device from tube. If flow is restored, check inhalation device. Replace it if necessary. 5 Unscrew humidifier (or outlet nipple). If flow is restored, check that humidifier is not obstructed. Check that the humidifier seal is correctly placed and in good condition. Correctly handtighten humidifier bottle.
• Autonomy of portable canister is less than usual	• Portable canister not completely full	Wait for hissing sound due to escaping gas to change before lifting vent lever.
• Liquid O ₂ leakage during disconnection portable - stationary	• Portable misconnected to the main canister • Blocks of ice gathered on the stationary canister filler opening	Immediately reconnect the 2 canisters then disconnect them. If leakage continues, reconnect them and contact your distributor
• Impossible disconnection portable - stationary	• stationary canister filler opening frozen	preventive : dry the filling outlets before each connection. Wait for a few minutes warming up the try again to disconnect the portable canister.

For any other problem, call your distributor.

6. Symbols - Abbreviations

	:	Type B device
	:	Attention : consult the accompanying documentation
	:	Audible alarm cutoff
	:	not toxic and not flammable gas
	:	Gas activating combustion
	:	Do not smoke in the room.
	:	Keep away from flame.
	:	Do not grease.
	:	Battery charge indicator
	:	Keep in vertical position

CE 0459 : Accordance with the directive 93/42/CEE (14th June 1993), about medical devices.

7. Standards

- EN 60 601-1 (Ed. 1995) : Electromedical devices safety
- EN 60 601-1-2 (Ed. 1993) : Electromagnetic compatibility
- Directive CEE 93/42 about medical devices
- ADR : Carrying on the road

8. Technical safeties

- Main and portable canister :
 - Automatic closing check valve (filling)
 - Hoses assembly automatic locking device
 - Ambient air valve (canister opening during their filling)
 - Safety relief valves (1,5 bar and 2 bars) (natural evaporation in case of not use)
 - Relief valve heater (ambient temperature operation)

9. Material in direct or indirect contact with the patient

- Canister :
 - Stainless steel
 - Aluminium
- Filters :
 - Sentered brass
- Miscellaneous :
 - Polycarbonate
 - PTFE
 - Silicon
 - Brass

10. Operating principle

Limits of environmental conditions		
Operating temperatures	°C	de 0 to + 40
Relative humidity (operating)		de 30% to 75%
Atmospheric pressure (operating)	hPa	de 700 to 1060
Storage temperatures	°C	de -40 to + 70
Relative humidity (storage)		de 10% to 100%
Atmospheric pressure (storage)	hPa	de 700 to 1060

Based on liquid oxygen regulated evaporation principle , the **FREELOX** system releases medical gaseous oxygen at ambient temperature under a constant adjustable flowrate, from liquide oxygen packaged in variable capacity isolated canisters.

		Stationary canister			Portable	
		44 l	32 l	20 l	1.2l	0.5l
Diameter	mm	365	365	365		
High	mm	880	740	600	365	295
Length	mm	/	/	/	195	195
Width	mm	/	/	/	130	130
Weight (empty)	kg	25.5	22	18.5	2.2	1.6
Weight (full)	kg	77	58.5	41	3.5	2.1
Capacity O2 liquid	l	44	32	20	1.15	0.45
Gas equivalent (15°C)	m³	33.6	24.4	15.3	0.83	0.34
Maximal operating pressure	bar	1.5	1.5	1.5	1.4	1.4
Normal operating pressure	bar	1.4	1.4	1.4	1.4	1.4
Evaporation	l/day	0.6	0.6	0.6	0.5	0.5
Adjustable flow rate (± 10 %)	l/min	0 to 7	0 to 7	0 to 7	0 to 7	0 to 7
Filling time (average)		7 min	5 min	4 min	60 sec	45 sec

Installation checklist / After maintenance checklist (before each patient installation)Yes No
done

- Verify the good general condition of the stationary canister. 
- Check that the humidifier outlet is not distorted. 
- Check that the flow rate selector knob rotates easily and that the indexing is aligned correctly for all settings. 
- Verify operation of the ambient valve. Assure that the portable canister is easily placed on the stationary canister. 
- Connect portable on the stationary canister to insure good connection and to check vertical alignment of the assembly (portable on stationary canister). 
- Check that the level indicator indicates a coherent value. 
- Check battery indicator. Change battery if necessary. 

Starting checklist (before each use)

- connect a patient circuit that is clean and in good general condition. Change it if necessary. 
- Fill the humidifier with water and re-fasten it tightly in order to avoid all leaks. 
- Check the good general condition of tube between cannula and humidifier (flexibility, not obstructed, connection without leak, ...). Replace it if necessary. 
- Check the level of water in the water collecting bottle. Empty it if necessary (clean and dry it regularly). 
- Adjust the oxygen flow rate to the prescribed value. 
- Check that the nasal cannula are not obstructed and that oxygen flows out (cf test with water). 
- In case of use with portable canister outside, check that the oxygen quantity is sufficient for the planned travel. 

Declaration for home - fire (car) insurance

Name :
Address :

Compagny :
Policy n° :

For the attention of the Director

Dear sir,

I wish to bring to your attention that my state of health means that I require prolonged oxygen therapy.

This treatment therefore includes the use and storage of a liquid oxygen container in my home.

(This may require me on occasion to transport a 0.5l or 1.2l container of liquid oxygen in my car. This use is exceptional and strictly therapeutic).

Please advise me if I need to take out additional cover for this, or if, like most other insurance companies, your company agrees to cover this risk without any surcharge.

Yours faithfully.

MAINTENANCE SHEET**1st year**

Performed on :

Nb of hours :

Name of Technician :

Signature
and
stamp**2nd year**

Performed on :

Nb of hours :

Name of Technician :

Signature
and
stamp**FREELOX n°** :**Put into service on** :**Maintenance provided by** :**Your distributor** :**address** :**Telephone** :**3rd year****Complete maintenance**

Performed on :

Nb of hours :

Name of Technician :

Signature
and
stamp**4th year**

Performed on :

Nb of hours :

Name of Technician :

Signature
and
stamp

Preventive maintenance of equipment must be performed according to the manufacturer's instructions given in the maintenance manual and any updates.

The operations must be performed by technicians who have received the appropriate training.

Use only original replacement parts.

Upon request, the manufacturer will provide electrical and pneumatic diagrams, a component list, calibration procedures, and all other necessary information for the use and maintenance of this device.

5th year**Complete maintenance ADR**

Performed on :

Nb of hours :

Name of Technician :

Signature
and
stamp**6th year**

Performed on :

Nb of hours :

Name of Technician :

Signature
and
stamp**7th year**

Performed on :

Nb of hours :

Name of Technician :

Signature
and
stamp**8th year**

Performed on :

Nb of hours :

Name of Technician :

Signature
and
stamp**Taema**

*AIR LIQUIDE Healthcare : a mission, an ethic, an ambition.
In more than 40 countries, the mission of the AIR LIQUIDE
Healthcare employees is to improve healthcare extended to
patients,*

*whether in hospitals or at home. Their ethic requires commitment to
the welfare of patients and support to all people who look after
them.*

*Their ambition is the development of the AIR LIQUIDE Group in the
Healthcare profession.*

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**Taema S.A.**