Vario

SUCTION PUMP

- EN Instructions for use
- DE Gebrauchsanweisung
- FR Mode d'emploi
- IT Istruzioni per l'uso
- NL Gebruiksaanwijzing
- SV Bruksanvisning
- FI Käyttöohjeet
- NO Bruksanvisning
- ES Instrucciones de uso
- PT Instruções de utilização

Vario 18 AC/DC

Table of contents

Congratulations	4
1 Warnings and safety instructions	5
Safety instructions	7
2 Power supply and battery operation	8
Battery low (only AC/DC versions)	8
Battery empty (only AC/DC versions)	8
3 Description	8
Introduction	8
Intended use	8
Indications	9
Contraindications	9
Intended user	9
Intended Patient Population	9
Important note	9
4 Overview	10
Definition of vacuum	10
Versions of the pump	10
Front of the pump	10
Bottom of device	11
Operating elements and battery status	11
Acoustic indicator	11
On/off switch not c/i versions	11
On/off switch c/i versions	11
5 Installation	12
Check initial delivery	12
6 Preparation for use	13
Checks before use	13
Assembly of the basic configuration	13
Assembly of the REUSABLE collection system	14
Assembly of the DISPOSABLE collection system	15
7 Operating instructions	16
Connect Vario to mains power	16
Check maximal vacuum for functional check	16
Changing vacuum level	17
Placing out of operation after use	17
8 Changing reusable jar	18
9 Changing disposable liner	19

10	Troubleshooting.....	20
	Motor not running.....	20
	Insufficient vacuum.....	20
	No LED lit (AC/DC version only).....	20
	LED lit up (AC/DC version only).....	20
	LED blinks slowly (AC/DC version only).....	20
	LED blinks rapidly (AC/DC version only).....	20
11	Replacing defective fuse.....	21
12	Setting interval times.....	22
13	Replacing the QuatroFlex.....	23
14	Replacing the rechargeable battery.....	24
15	Battery data and test.....	25
	Duration of pumping.....	25
	Charging time (charged up to 100%).....	25
	Testing the rechargeable battery.....	25
16	Battery calibration.....	26
	Calibrating NiMH rechargeable battery.....	26
	Storage of the rechargeable battery.....	26
17	Filter test.....	27
18	Cleaning guidelines.....	28
	General notes.....	28
	Medela recommended surface cleaning agents for pump housing.....	28
	Water.....	28
	Cleaning/ disinfection machines.....	28
	Disposable products.....	28
19	Accessories overview.....	30
20	Accessories legend.....	31
21	Technical specifications.....	33
22	Signs and symbols.....	34
23	Warranty and maintenance.....	36
	Warranty.....	36
	Maintenance/routine check.....	36
24	Disposal.....	36
25	Technical documentation.....	37
	Technical documentation/ EMC.....	366

Congratulations

By choosing the Vario, you have acquired a high-quality suction pump which is adaptable to your needs. The QuatroFlex suction aggregate is the innovative continuation of the proven Medela piston/cylinder system. As all Medela suction pumps, this pump provides reliable and simple suction. Its simple handling and cleaning as well as the safety features are additional advantages you receive. A comprehensive range of accessories makes the Vario ideally suited to a wide range of medical applications and can be used for continuous operation. Contact us – we will be pleased to advise you.

1 Warnings and safety instructions

WARNINGS

Indicates a potentially hazardous situation which, if not avoided, could result in death or serious injury.

CAUTIONS

Indicates a potentially hazardous situation which, if not avoided, could result in minor or moderate injury.

Safety related tip

Indicating useful information about the safe use of the device.

The Vario is approved exclusively for the use as described in these instructions for use. Medela can only guarantee the safe functioning of the system when the Vario is used in combination with the original Medela accessories (collection system, tubings, filters etc. – see chapters 19/20).

The Vario is EMC-tested in conformity with the requirements of IEC 60601-1-2:2007 and can be used in the vicinity of other EMC-tested devices that fulfil the requirements as outlined in the IEC 60601-1-2 standard. Untested HF sources, radio networks or the like can impair the function of the device and should not be operated in combination with the Vario.

Please read and observe these warning and safety instructions before operation. These instructions for use must be kept with the device for later reference.

Please note that these instructions for use are a general guide for the use of the product. Medical matters must be addressed by a physician.

Medela considers herself only responsible for the effect on BASIC SAFETY, reliability and performance of the Vario if it is used in accordance with the instructions for use.

Caution

Federal US law restricts this device to sale by or on the order of a health professional.

Subject to change.

WARNINGS

- For use only by medically trained persons who have been adequately trained in suction procedures and in the use of aspirators.
- For patients for whom a breakdown of the device can lead to a critical situation (e.g. patients with acute dyspnoea or severe catarrhal congestion), a replacement device must always be available.
- The device must not be used for suctioning explosive, easily flammable or corrosive liquids
- The connecting tubing supplied with the device must never come into direct contact with the suction area. A sterile suction catheter must always be used (risk of infection).
- Before cleaning the device, pull the plug out of the fixed mains socket.
- No modification of this equipment is allowed.
- Do not connect this device to a passive drainage tube.
- Consult the indications for use and consider risk factors and contraindications before using the Vario. Failure to read and follow all instructions in this manual prior to use may result in serious or fatal injury of the patient.

CAUTIONS

- Incorrect use can cause pain and injury to the patient.
- Do not use sterile accessories when the sterile packaging is damaged.
- The use of mobile telephones, LAN/WLAN, walkie-talkies (two-way radios) and cordless telephones sets can affect the Vario pump. A safety distance of min. 3.3ft (1 m) to the Vario pump is recommended.
- Portable and mobile RF communications equipment can affect medical devices.
- The patient should be monitored regularly according to the physicians instructions and facility guidelines. Objective indications or signs of a possible infection or complication must be met immediately (e.g. fever, pain, redness, increased warmth, swelling or purulent discharge). Non-observance can lead to considerable danger of the patient. Monitor the Vario frequently for operating status.

Safety related tip

- The Vario suction pump is Magnetic Resonance (MR) unsafe. Do not take the pump into the MR environment.
- For safety tests, it is assumed that the device is serviced and repaired throughout its service life in accordance with the service manual. Because the Vario is a device of safety class II (EN IEC 60601-1, A1, A2:1995), the safety tests are confined to visual inspection of the housing and power cord. This test must be carried out before each use.
- Testing the patient leakage current: Please see service manual for details.
- The protection of the Vario against the effects of the discharge of a cardiac defibrillator is dependent upon the use of appropriate cables.
- Separation from the mains is only assured through the disconnection of the mains plug and the fixed mains socket.
- Before you plug in the device, please check that your local power supply is the same as the voltage given on the specification plate.

Safety instructions

- The Medela Vario is a medical device that requires special safety measures in regard to EMC. It must be installed and put into operation in accordance with the attached EMC information.
- In the case of overflow, inform the internal technical service immediately and perform the tasks in the service manual.
- In each of the following cases, the device must not be used and it must be repaired by Customer Services:
 - if the power cord or the plug are damaged
 - if the device is not functioning perfectly
 - if the device is damaged
 - if the device shows clear safety defects.
- Keep the power supply cord away from hot surfaces.
- The mains plug and the on/off switch must not come into contact with moisture.
- Never pull the mains plug out of the fixed mains socket by pulling on the power supply cord!
- Never leave the device unattended when it is switched on.
- The pump must stand upright during use.
- Never place the power supply cord around your neck.
- Never use the device at high room temperatures, while bathing or showering, if you are very tired or in an environment where there is a risk of explosion.
- Never place the device in water or other liquids.
- When using single use, sterile products, please note that they are not intended to be reprocessed. Reprocessing could cause loss of mechanical, chemical and/or biological characteristics. Reuse could cause cross contamination.
- Contact your local Medela customer service representative for assistance with product operations.
- Do not place a 2l suction jar in the lateral adaptor.
- The overflow protection/bacteria filter protects the pump against overflow. Caution: Suctioning is interrupted when the filter is wet/moist or clogged. Test the filter before each use and during use periodically (see chapter 17).

These instructions for use must be kept for later reference.

2 Power supply and battery operation

The Vario is a mains-powered suction pump. Before you plug in the device, please check that your local power supply is the same as the voltage given on the specification plate. If the Vario AC/DC versions are running continuously, a maximum of 30 minutes of battery time is expected before it is necessary to recharge. If the Vario AC/DC version is not in use, the battery must be calibrated every 60–90 days.

Battery low (only AC/DC versions)

During battery operation a slow blinking of the LED indicator and an acoustic signal (beep - beep - beep) indicates, that the rechargeable battery is soon discharged. The function of the Vario is maintained for 10 minutes (20 minutes Vario low vacuum versions), but the battery must be recharged as soon as possible.

Battery empty (only AC/DC versions)

As soon as the rechargeable battery is empty, the pump comes to a stop. The proper function of the pump is no longer maintained. The rechargeable battery must be recharged.

3 Description

Introduction

The Vario is a high-quality suction pump. It is powered by the well-proven QuatroFlex system and guarantees maximum suction performance for many suctioning needs. It ideally combines easy handling and cleaning with safety features to ensure optimal operation. You can choose from a comprehensive range of accessories from Medela to configure the pump to many medical applications. The device should not be used adjacent to or stacked with other equipment. If adjacent or stacked use is necessary, the device should be observed to verify normal operation in the configuration in which it will be used.

Intended use

The Vario 8/18/ci suction pumps are indicated for aspiration and removal of surgical fluids, tissue (including bone), gases, bodily fluids (including vomit) or infectious materials from a patient's airway or respiratory support system, either during surgery or at the patient's bedside. Generally the Vario 8/18/ci is intended to be used for a variety of suctioning procedures including nasopharyngeal, tracheal, surgical, gastrointestinal suctioning either in «constant» or «intermittent» mode. Please make sure you have a correct Vario version for the required application.

Indications

Gastric suction: Only with Vario 18 c/i: When the device is used for gastric suction, the patient's condition and correct functioning of the pump (on the basis of acoustic and visual signs) must be checked at least once every four hours.

ENT indications: For patients for whom a breakdown of the device can lead to a critical situation (e.g. patients with acute dyspnoea or severe catarrhal congestion), a replacement device must always be available.

Contraindications

Not suitable for setting at a low vacuum, as needed for example for thoracic drainage without specialized accessories.

Intended user

The Vario should only be operated by properly trained staff. These persons must not be hard of hearing or deaf and must have adequate visual faculty. The training has to be refreshed at least once a year.

Intended Patient Population

The Vario is intended to be used on patients only exhibiting conditions as described in the indications for use.

Important note

Compliance with proper surgical procedures and techniques is the responsibility of the physician. Each physician must evaluate the appropriateness of the treatment based on his own knowledge and experience.

Service life

The service life of the device is three years.

4 Overview

Definition of vacuum

By the application of medical aspiration devices, vacuum is normally given as the difference (in absolute figures) between absolute pressure and atmospheric pressure or as negative values in Kilopascal (kPa). In this document, the indication of -10 kPa for example always refers to a pressure range in kPa below atmospheric ambient pressure (according to EN ISO 10079:1999).

Versions of the pump

Vario 8 AC

Vario 8 AC/DC

Vario 8 c/i AC/DC

Vario 18 AC

Vario 18 AC/DC


Vario 18 c/i AC/DC


Front of the pump

Bottom of device

Operating elements and battery status

not lit Battery is fully charged

lit up – Battery is being charged
– Pump is connected to the mains


blinks slowly Battery is soon empty

blinks rapidly Battery is defective


Acoustic indicator

During battery operation an acoustic signal (beep - beep - beep) indicates that the rechargeable battery is soon discharged.

On/off switch not c/i versions

On/off switch c/i versions

5 Installation

Check initial delivery

Check the delivery package of the Vario for completeness and general condition.

Vario portable version

REF	Vario 8 AC	026.XXXX
	Vario 8 AC/DC	026.XXXX
	Vario 8 c/i AC/DC	026.XXXX
	Vario 18 AC	026.XXXX
	Vario 18 AC/DC	026.XXXX
	Vario 18 c/i AC/DC	026.XXXX

Silicone Tubing $\varnothing 7 \times 12$ mm, 25 cm with 2 coupling pieces

REF 077.0931

5 x Disposable overflow protection/bacteria filter with Medela connections

REF 077.0571

Instructions for use

REF 177.2604

6 Preparation for use

WARNINGS

For use only by medically trained persons who have been adequately trained in suction procedures and in the use of aspirators.

CAUTIONS


- The Vario must remain in an upright position during use.
- Sterile accessories must be checked on the integrity of the packaging before use.
- Non sterile and reusable accessories must be cleaned, disinfected and/or sterilized according to the Medela cleaning guide (product code 200.2391).

6.1 Checks before use

- Check the Vario system before use for damage of the power cord or plug, obvious device damage or safety defects and proper functioning of the device.
- Check for completeness and general condition of the Vario delivery package.
- Only for AC/DC versions: make sure that the rechargeable battery is calibrated according to chapter 16.
- Check all accessories prior to use:
 - suction jars, lids and liners for cracks, brittle and flawed spots. Replace if necessary.
 - Tubing for cracks, brittle areas and that connectors are firmly attached. Replace if necessary.
 - As an additional safety test, evacuate the system (including jars) to maximum vacuum several times before actual use.

6.2 Assembly of the basic configuration


1. 
 - 1.1 Attach the overflow protection bacteria filter to the Vario with the arrow pointing in the flow direction. The filter must be exchanged after an overflow and after each patient, but latest after a week.
2. Attach all necessary accessories according to your needs. See the system overview for more details.


6.3 Assembly of the REUSABLE collection system


1. If you use a DISPOSABLE collection system proceed with step 6.4.

2.1 Attach mechanical overflow protection to lid.
2.2 Pull gently downwards to make sure it is open/deactivated.

3.1 Attach the lid to the jar.
3.2 Lock with the two lid clamps.

4.1 Attach the jar to the Vario.
4.2 Connect tubing to the filter.
4.3 Connect tubing to the lid of the jar (marked with «vacuum»).

5.1 Connect patient tubing to the lid of the jar (marked with «patient»).

6.4 Assembly of the DISPOSABLE collection system


positions for jars:

1. If you use a REUSABLE collection system proceed with step 6.3.

2.  2.1 Attach the Vario jar adapter to the Vario.

3.  3.1 Choose the 1.5l size (yellow code) for the Vario.
3.2 Prepare and insert corresponding liner size (yellow) according to instructions delivered with liners. The disposable liner has an integrated overflow protection filter. No action is necessary to activate it.

4.  4.1 Attach the jar to the Vario.
4.2 Connect tubing to the filter.
4.3 Connect tubing to the jar. The tubing connector is on top of the jar.

5.  5.1 Connect patient tubing to the lid of the liner using either the coloured angle piece or the patient port directly (depending on tubing size).

7 Operating instructions

CAUTIONS

The Vario is to be set up in such a way, that a separation from the mains supply can be easily managed.

7.1 Connect Vario to mains power

1. Check the pump before use following the instruction in chapter 6.1.

2. **All versions: if the pump is connected to a fixed mains socket:**

Plug in the mains plug of the power cord to a fixed mains socket.

or

- AC/DC versions: if the pump is operated with a 12VDC power source:**


Connect the optional 12VDC connection cable to the socket on the side of the pump and connect it to the 12VDC power source.

or


- AC/DC versions: if the pump is operated with the rechargeable battery:**

No connection to a power source is needed. Make sure that the battery was charged for five hours.

7.2 Check maximal vacuum for functional check

1. Switch on the Vario. The pump will start running. Caution: **c/i versions** must be switched on in «constant» mode.

- 2.1 Press and turn the vacuum regulator to the right to set maximum vacuum.

- 3.1 Seal the end of the patient tubing with your thumb.
3.2 Compare the maximum vacuum according to the specification. See chapter 10 if the maximum vacuum is not reached.

Specifications:

Altitude above sea level:	Vario 8 / Vario 8 c/i Max. Vacuum: low vacuum	Vario 18 c/i Max. Vacuum: medium vacuum	Vario 18 Max. Vacuum: high vacuum
▼	▼	▼	▼
+ 2000m	- 7.0 kPa	- 43 kPa	- 59 kPa
+ 1000m	- 7.9 kPa	- 48 kPa	- 66 kPa
+ 500m	- 8.5 kPa	- 51 kPa	- 70 kPa
0m	- 9.0 kPa	- 55 kPa	- 75 kPa


(Tolerance +/- 10 %)

7.3 Changing vacuum level

1. Clamp patient tubing
- 1.2 Push and turn vacuum regulator to select the correct vacuum according to the particular application.
- 1.3 Check vacuum gauge for setting.

7.4 Placing out of operation after use


1. **1.1 Switch the on/off switch of the Vario to position «0» (c/i - versions: «Off»).**
2. **If the pump is connected to a fixed mains socket:** Disconnect the mains plug from the fixed mains socket.
or
If the pump is connected to a 12VDC power source: Disconnect the 12VDC plug from the 12VDC power source.
3. Clean and disinfect the Vario. See chapter 18.

8 Changing reusable jar

1.

- 1.1 Switch the on/off switch of the Vario to position «0» (c/i - versions: «Off»).
- 1.2 Remove patient and vacuum tubings from the lid.

2.

- 2.1 Remove full jar from the Vario.

3.

- 3.1 Prepare new jar and unpack it.
- 3.2 Secure new jar to the Vario
- 3.3 Reconnect vacuum tubing and new patient tubing firmly.

4. Empty jar and dispose of patient tubing in accordance with local guidelines and inhouse guidelines.
5. Switch on the Vario. The pump will start running.
6. Set vacuum according to the particular application. See chapter 7.3.

9 Changing disposable liner


WARNINGS

Disposable liners are not intended to be reprocessed. Reprocessing could cause loss of mechanical, chemical and/or biological characteristics.


1.5l

1.

1.1 Remove patient tubing with colored angle piece from lid.

1.2 Plug patient port.

2.

2.1 Switch the on/off switch of the Vario to position «0» (c/i - versions: «Off»). Reduce vacuum to minimum by pressing and turning vacuum regulator to the left.


3.

3.1 Remove full liner from the jar.

4. Dispose of the liner and patient tubing in accordance with local guidelines and inhouse guidelines.

5. Prepare new liner, unpack it and insert it in the jar. Press down lid.

6.

6.1 Connect new patient tubing firmly.

7. Switch on the Vario. The pump will start running.

8. Set vacuum according to the particular application. See chapter 7.3.

10 Troubleshooting

Motor not running

Check if:

- the Vario is switched on.
- the mains voltage is correct and if the mains plug is inserted correctly into the fixed mains socket.
- the internal battery is charged (only AC/DC version under battery operation)
- the fuse on the bottom of the Vario is not defective. For replacing the defective fuse see chapter 11.

Insufficient vacuum

Check if:

- the vacuum regulator is set correctly.
- the tubings are not defective or broken. If necessary, replace.
- all plug-in connections are tight.
- the overflow protection is deactivated/open. If the overflow protection is activated, deactivate as shown under 6.3a/2.2.
- the suction jar/lid have no cracks, brittle areas, discolouration. Replace if necessary.
- the disposable system has no cracks, brittle areas, discolouration. Replace if necessary.
- filter is not clogged. To test if the filter is clogged see chapter 17.
- the QuatroFlex is defective. To replace the QuatroFlex see chapter 13.

No LED lit (AC/DC version only)

The rechargeable battery is fully charged.

LED lit up (AC/DC version only)

- The rechargeable battery is being charged.
- The pump is connected to the mains.

LED blinks slowly (AC/DC version only)

The rechargeable battery is soon empty.

LED blinks rapidly (AC/DC version only)

The rechargeable battery is defective. For replacing the defective battery see chapter 14.

11 Replacing defective fuse

WARNINGS

Before replacing the fuse, pull the mains plug from Vario out of the fixed mains socket.

1. Disconnect the Vario from the mains.

- 2.

- 2.1 Open (turn left) the fuse holder on the bottom of the Vario.

3. Replace the defective fuse(s). Make sure that the technical specification are corresponding.

AC-versions: **T 0.8 AL, 5x20, 230–240V**

T 1.25 AL, 6.3x20, 120V

AC/DC-versions: **T 1.0 AL, 5x20, 100–240V**

- 4.

- 4.1 Close (push back and turn right) the fuse holder.

5. Reconnect the Vario to the mains again.

6. Switch on the Vario again.

12 Setting interval times


(c/i versions only)

WARNINGS

Before setting the interval times, switch off the Vario and pull the mains plug from Vario out of the fixed mains socket.


1.

1.1 Unscrew the 4 screws on the back of the Vario.

1.2 Open the Vario by removing the back cover.

2.

2.1 Determine the location of the red jumper.

2.2 Place the red jumper according to the desired interval:

On-time	Off-time
8 sec.	4 sec.
16 sec.	8 sec.
32 sec.	16 sec.

3. Reattach the back cover to the Vario and tighten the 4 screws.





4. Plug the Vario to a fixed mains socket and test the correct functioning of the Vario (see chapter 7).

13 Replacing the QuatroFlex

WARNINGS

Before replacing the QuatroFlex, switch off the Vario and pull the mains plug from Vario out of the fixed mains socket.

1. 
 - 1.1 Unscrew the 4 screws on the back of the Vario.
 - 1.2 Open the Vario by removing the back cover.
2. 
 - 2.1 Remove the 3 marked tubings from the QuatroFlex.
3. 
 - 3.1 **Turn** the QuatroFlex 45° to the left.
 - 3.2 **Pull out** the QuatroFlex and remove it from the Vario.
4. Clean inside of the Vario.
5. 
 - 5.1 **Insert** new QuatroFlex.
 - 5.2 **Turn** it 45° to the right.
 - 5.3 **Reattach** the 3 tubings to the QuatroFlex.
6. Reattach the back cover to the Vario and tighten the 4 screws.
7. Plug the Vario to a fixed mains socket and test the correct functioning of the Vario (see chapter 7).

14 Replacing the rechargeable battery


(AC/DC versions only)

WARNINGS

Before replacing the rechargeable battery, switch off the Vario and pull the mains plug from Vario out of the fixed mains socket.

1.

1.1 Unscrew the 4 screws on the back of the Vario.


1.2 Open the Vario by removing the back cover.

2.

2.1 Detach the battery cable by pulling the plug.

3.

3.1 Remove the old rechargeable battery.

3.2 Insert the new rechargeable battery.

4.

4.1 Reattach the battery cable.

5. Reattach the back cover to the Vario and tighten the 4 screws.

6. Plug the Vario to a fixed mains socket and calibrate the rechargeable battery (see chapter 16).

15 Battery data and test

(AC/DC versions only)

CAUTIONS

Only charge the rechargeable battery when the acoustic signal sounds and the charging indicator blinks slowly. The rechargeable battery is charged during mains operation.

Duration of pumping

If the battery is fully charged, the operating time is about:


- low vacuum versions: approx. 120 minutes
- high and medium vacuum versions: approx. 30 minutes.

Charging time (charged up to 100%)


- low vacuum versions: approx. 5 hours
- high and medium vacuum versions: approx. 5 hours.

Testing the rechargeable battery

1. Make sure that the rechargeable battery is fully charged.
2. Pull out the mains plug from the fixed mains socket.

- 3.1 Press vacuum regulator and turn it to the right to select maximum vacuum.

4. Switch on the Vario. The pump will start running. Caution: **c/i versions** must be switched on in «constant» mode.
5. Measure the time until the pump comes to a stop. The battery is ok, when the following values are reached:
 - low vacuum versions: > 120 min.
 - high and medium vacuum versions: > 30 min.

If the values are not reached, repeat test or calibrate the battery (see next page).

16 Battery calibration

(AC/DC versions only)

CAUTIONS


Calibrate rechargeable batteries before the first use, or after a storage period >2 months. Replace the battery if the calibration process fails.

Calibrating NiMH rechargeable battery

1. Make sure that the battery is fully charged.

2. Pull out the mains plug from the fixed mains socket.

3.

3.1 Switch on the Vario. The pump will start running. Caution: **c/i versions** must be switched on in «constant» mode.

4. Let the pump run until the rechargeable battery is empty and the pump comes to a stop.

5. Plug the mains cable into a fixed socket. The rechargeable battery is being charged.




6. Wait until the charging indicator goes out.

7. Repeat step 1 to 6. The Vario AC/DC is then ready to use.

Storage of the rechargeable battery

In order to keep the self-discharge at a minimum, store the Vario and the replacement rechargeable batteries at temperatures below 25°C (77 °F). Repeat the calibrating process every 60–90 days.

17 Filter test

1.  1.1 Remove tubing from filter.
2.  2.1 Select maximum vacuum.
3.  3.1 Switch on the Vario. The pump will start running. Caution: **c/i versions** must be switched on in «constant» mode.
4.  4.1 Read the vacuum.
4.2 Replace the filter if the vacuum exceeds the values given in the following table:


Vario 8 / Vario 8 c/i low vacuum	– 3kPa
Vario 18 high vacuum	– 10kPa
Vario 18c/i medium vacuum	– 10kPa

18 Cleaning guidelines

WARNINGS

After each use, the parts that had been in contact with the aspirated secretions are to be cleaned, disinfected, sterilised or disposed of **according to the table on the next page**.

Before cleaning the device, pull the mains plug out of the fixed mains socket.

General notes

- These are general recommendations only that may be adjusted individually, based on the hospital's specific directives and cleaning practices and policies.
- Detailed information according to Medela cleaning instructions (product code 200.2391).
- Wear protective gloves for cleaning/disinfection.
- Dispose of fluids such as blood and secretions and the parts contaminated with them in accordance with internal hospital guidelines.

Medela recommended surface cleaning agents for pump housing

- **Hexaquart plus** (aldehyde-free), B. Braun Medical AG, www.bbraun.ch
- **Perform** [concentrate], Schülke & Mayr, www.schuelkemayr.com

Water

Use only the purest quality of water for cleaning. Water hardness is a serious consideration since deposits left on medical products may not be properly decontaminated. Use deionised water in order to reduce this problem. The final rinse water is bacterial free and contains no endotoxins.

Cleaning/disinfection machines

Can be used to desinfect parts from the table on the next page. A hot water rinse (maximum temperature 100 °C) may provide a medium-to high level of disinfection. Every section of the constituent parts must be accessible in order to ensure efficient cleaning. We recommend using a cleaning/disinfecting machine that has been approved by the Robert Koch- Institute and complies with ISO 15883. Recommended temperature for noncritical medical devices (i.e. those that only come into contact with uninjured skin) is 90 °C for 1 minute. The time is increased to 5 minutes for all medical devices that are considered to be critical.

Disposable products

These are single use products not intended to be reprocessed. Reprocessing could cause loss of mechanical, chemical and/or biological characteristics. Reuse could cause cross contamination.


PSU* jars, lids	PC** jars	Silicone tubing	Pump housing Mains cable	Plastic accessories	Legend: * Polysulfone ** Polycarbonate
x	x	x	x	x	1. Disassembly Separate all individual parts before cleaning, disinfecting and sterilising.
x	x	x			2. Cleaning Clean components in hot water (60–70 °C) containing a detergent with a pH range between 6.0 and 8.0 only, in order to avoid damaging the instruments and containment devices. Enzymatic detergents help with the removal of organic matter, such as blood. Detergents should be used according to their manufacturer's recommended concentration levels. Some alkaline detergents have been formulated to be safe for reprocessing medical devices. The manufacturers of these detergents should provide information about specific materials that may be damaged by their detergents. Soak all parts thoroughly with warm, soapy water or in enzymatic detergent for 1–5 minutes. 1. Remove visible dirt with a cleaning tool – general purpose cleaning brushes, such as pipe cleaners or non-abrasive lint cloths. Brushes and pipe cleaners should fit snugly but still be able to be moved around easily within the area to be cleaned. Rinse thoroughly in clear water. Allow to dry 2. Check the instruments for visible dirt and repeat these steps if necessary.
			x	x	3. Cleaning Wipe with detergent mentioned above. Use clean cloth to dry.
x	x	x		x	4. Disinfect Soak at room temperature for 30 minutes in a disinfection/ cleaning solution (A). After the reaction time, use water to rinse the solution residue from the individual parts and rinse the tubing. Rinse for at least 2 minutes with the cleaning solution (B). (A) fluid aldehyde-free disinfection solution with cleaning effect for use as a soaking bath, bactericide, fungicide, limited virucide, with good material compatibility for stainless steel, non-ferrous material and plastics including silicone, slightly alkalic. e.g. 2 % neodisher® Septo MEDsolution (B) Fluid, pH-neutral, enzymatic cleaning solution for the treatment of instruments by machine or manually with very good material compatibility for stainless steel, non-ferrous material and plastics including silicone. or use cleaning/ disinfecting machine. or use your dishwasher at home.
x		x			5. Sterilise Remove any chemical disinfectant residue before autoclaving. Do not stack during autoclaving. In the autoclave at 134 °C for 10 minutes. The sterilisation device must comply to ISO 17665:2003. The following procedures are acceptable: In the autoclave with saturated steam at 132 °C, with triple pre-evacuation, sterilisation time of 6 minutes, or at 135–137 °C in an porous load vacuum autoclave for 3–3.5 minutes
x	x	x		x	6. Storing components After disinfection/ sterilising, store the components in sterile foil until required for use.

19 Accessories overview

WARNINGS

Vario was verified in combination with the accessories listed below. For a correct and safe operation use Vario with these accessories only. Further information is supplied on the instruction sheet of the individual accessory.

Patient tubings/tubing connectors		<i>without coupling pieces, 7 x 12 mm</i>	<i>with 1 coupling piece, 7 x 12 mm</i>
 077.1020	 077.1017	 077.0900: 25 m 077.0901: 1 m 077.0902: 2 m 077.0905: 5 m	 077.0912: 0.6 m 077.0913: 1 m
Liners		Reusable lids	
 077.0083/84	 077.0440	 077.0420/30	Fingertip tubing  077.0951
PC jars		Reusable jars	
 077.0082	 077.0105	 077.0125: 0.25 l 077.0155: 0.5 l 077.0110: 1 l	
Connectors and Tubings (max. length 60 cm)		Reusable sets	
 077.1019	 077.1018	 077.1015/16	 077.0715: 1 l
Filters			
 077.0092	 077.0571	 077.0577	
* Drainage valve	* Specimen cups		
 077.1022	 077.0562/64 077.0094		
Additional accessories			
 077.0823	 077.0825	 077.0821	 077.0831
 077.0832			
			
Vario 8 AC 026.XXXX Vario 8 AC/DC 026.XXXX Vario 8 c/i AC/DC 026.XXXX Vario 18 AC 026.XXXX Vario 18 AC/DC 026.XXXX Vario 18 c/i AC/DC 026.XXXX			

**Safety related tip**

If the pump is used together with non-Medela applied parts, they must:

- have CE mark and (if necessary) local registration
- be able to be attached to Medela accessories safely without impacting the performance of the pump.

Note: when combining Medela parts and non-Medela parts, you take over the responsibility of the entire system.

20 Accessories legend

REF

026.XXXX Vario 8 AC
 026.XXXX Vario 8 AC/DC
 026.XXXX Vario 8 c/i AC/DC
 026.XXXX Vario 18 AC
 026.XXXX Vario 18 AC/DC
 026.XXXX Vario 18 c/i AC/DC

Filters

077.0092 Disposable overflow protection/bacteria/odour filter with Medela connections
 077.0571 Disposable overflow protection/bacteria filter with Medela connections
 077.0577 Disposable overflow protection/bacteria filter with Medela and conical connections

Reusable jars

077.0110 Suction jar, polysulfone, 1 l
 077.0125 Suction jar, polysulfone, 0.25 l
 077.0155 Suction jar, polysulfone, 0.5 l

Drainage valve

077.1022 Drainage valve (for vacuum ports of reusable jars)

Liners

077.0083 Disposable suction liner 1.5 l
 077.0084 Disposable suction liner 1.5 l with solidifier

Reusable sets

077.0715 Reusable 1 l set with PSU suction jar

PC jars

077.0082 Suction jar PC, 1.5 l
 077.0105 Vario jar adapter

Specimen cups	077.0094	Specimen cup for disposable system. Can be attached to lid of liner (patient port) of Disposable Collection System
	077.0562	Disposable specimen cup, Ø6–10 mm
	077.0564	Disposable specimen cup, Ø 10–14 mm. Attach to lid of Reusable Collection System (patient port)
Reusable lids	077.0440	Small lid with conical patient connection, Ø6–10 mm and overflow protection device
	077.0420	Large lid with conical patient connection Ø6–10 mm and overflow protection device
	077.0430	Large lid with conical patient connection Ø 10–14 mm and overflow protection device
Patient tubings/tubing connectors	077.0902	Silicone tubing Ø7 x 12 mm, w/o coupling pieces, 2 m
	077.0912	Silicone tubing Ø7 x 12 mm with 1 coupling piece, 60 cm
	077.0900	Silicone tubing Ø7 x 12 mm, w/o coupling pieces, 25 m
	077.0901	Silicone tubing Ø7 x 12 mm, w/o coupling pieces 1 m
	077.0905	Silicone tubing Ø7 x 12 mm, w/o coupling pieces, 5 m
	077.0913	Silicone tubing Ø7 x 12 mm with 1 coupling piece, 100 cm
	077.0951	Disposable (PVC) tubing, 180 cm, with fingertip, sterile (applied part)
	077.1017	Double conical coupling piece. For connecting 2 tubings to each other
	077.1020	Y-Piece. For connecting 3 tubings to each other
	077.0921	Silicone tubing Ø7 x 12 mm with 2 coupling pieces, 30 cm
Additional accessories	077.0922	Silicone tubing Ø7 x 12 mm with 2 coupling pieces, 60 cm
	077.0931	Silicone tubing Ø7 x 12 mm, 25 cm, with 1 coupling piece and 1 coupling piece 90°
	077.0821	Rail holder
	077.0823	Universal holder
	077.0825	Rail holder
	077.0831	Carrying bag
	077.0832	Car connection cable for 12VDC

21 Technical specifications

Vario 8 / Vario 8 c/i:

low vacuum, - 9 kPa/- 68 mmHg (Tolerance: +/- 10 %)

Vario 18 c/i:

medium vacuum, - 55 kPa/- 413 mmHg (Tolerance: +/- 10 %)

Vario 18:

high vacuum, - 75 kPa/- 563 mmHg (Tolerance: - 10 %)

Measured at 0 m, atmospheric pressure: 1013.25 hPa Please note: vacuum levels may vary depending on location (meters above sea level, atmospheric pressure and temperature).

Vario 8: 8 l/min. (+/- 10 %)

Vario 18: 18 l/min. (+/- 10 %)

AC 3.5 kg

AC/DC 4.2 kg

Without jar

230–240 V, 50 Hz, 90 VA
230–240 V, 60 Hz, 90 VA
120 V, 60 Hz, 70 VA

100–240 V
50/60 Hz
80 VA

ISO 9001

ISO 13485

CE (93/42/EEC), IIa

h x w x d

380 x 170 x 285 mm

Transport/Storage

Transport / Storage

Rechargeable battery

Operation

IP21

22 Signs and symbols

0123 This symbol indicates the compliance with the essential requirements of the Council Directive 93/42/EEC of 14 June 1993 concerning medical devices.

This symbol indicates the class of the pump.

This symbol indicates a class II device.

IP21

This Symbol indicates the protection against ingress of solid foreign objects and against harmful effects due to the ingress of water.

This symbol indicates the compliance with additional USA and Canada safety requirements for medical electrical equipment.

This symbol indicates a type CF applied part.

This symbol indicates the legal specifications of the pump.

This symbol indicates MR UNSAFE.

This symbol indicates to follow instructions for use.

This symbol indicates that interferences may occur in the vicinity of equipment marked with this symbol.

This symbol indicates to consult instructions for use.

This symbol indicates the manufacturer.

This symbol indicates a CAUTION or WARNING associated with the device.

This symbol indicates the date of manufacture (four digits for the year and two digits for the month).

This symbol indicates safety related tip.

This symbol indicates that the device should not be used after the end of the year and month shown.

This symbol indicates a single use device. Do not reuse the device.

This symbol indicates manufacturer's catalogue number.

This symbol indicates manufacturer's serial number.

This symbol indicates manufacturer's batch code.

This symbol indicates the device is sterilized using ethylene oxide.

This symbol indicates the temperature limitation for operation, transport and storage.

This symbol indicates the humidity limitation for operation, transport and storage.

This symbol indicates the atmospheric pressure limitation for operation, transport and storage.

This symbol indicates to keep the device dry.

R_x only

This Symbol indicates a Prescription Device. Federal Law restricts this device to sale by or on the order of a licensed healthcare practitioner. (for US only).

This symbol indicates do not use the device if package is damaged.

This symbol indicates the maximum vacuum level of the pump.

This symbol indicates the flow levels of the pump.

pcs

This symbol indicates number of items.

This symbol indicates the number of items n that the content is sufficient for.

This symbol indicates the electrical specifications of the pump.

This symbol indicates the fuse

This symbol indicates do not dispose the device together with unsorted municipal waste (for EU only).

This symbol indicates alternating current.

This symbol indicates that the material is part of a recovery/recycling process.

This symbol indicates direct current.

This symbol indicates a carton package.

This symbol indicates the weight of the pump.

This symbol indicates to keep the device away from sunlight.

This symbol indicates the dimensions (h x w x d) of the pump.

This symbol indicates to handle the fragile device with care.

23 Warranty and maintenance

Warranty

Medela AG warrants the device will be free from defects in materials and workmanship for a period of 2 years (6 month for the rechargeable battery and the QuatroFlex) from the date of delivery. Faulty material will be replaced free of charge during this period if not resulting from abuse or misapplication. This will not apply to parts subject to wear and tear in use. To ensure compliance with this warranty as well as optimum service from Medela products, we recommend the exclusive use of Medela accessories with our appliances.

In no event shall Medela AG be liable for claims which exceed the scope of warranty described including liability for consequential damages, etc. The right to the replacement of faulty parts will not be recognized by Medela if any work has been carried out on the pump by unauthorized persons. This warranty is subject to the appliance being returned to a Medela service centre.

Maintenance/routine check

Routine checks and service work are only to be carried out by positions authorised by Medela. The routine check is to be carried out 1 x per year. (see service manual). The service manual, wiring diagrams, and descriptions are available upon request from Medela.

24 Disposal

The Vario comprises metals and plastics and should be disposed of in accordance with the European directives 2002/95/EC and 2002/96/EC. Additional, local guidelines must also be observed. With the AC/DC versions, the electronic components and the rechargeable battery must be disposed of separately, in accordance with the local regulations.

Please take care that you dispose the Vario and its accessories in accordance with your local disposal guidelines.

User information for the disposal of electrical and electronic equipment

This symbol means that the electrical and electronic equipment must not be disposed as normal household refuse. A correct disposal of this device protects and prevents possible damage to the environment or human health. For more information about the disposal contact the manufacturer, your local caregiver or healthcare provider. This symbol is only valid in the European Union. Please respect the relevant state laws and rules in your country for the disposal of electrical and electronic equipment.

25 Technical documentation

WARNINGS

Do not use other accessories than those specified or sold by the manufacturer as replacement parts for internal components as it may result in increased emissions or decreased immunity of the Vario pump.

EMC

Vario is EMC-tested in conformity with the requirements of IEC 60601-1-2:2007 and can be used in the vicinity of other EMC-tested devices that fulfill the requirements of the relevant IEC 60601-1-2 standard. Untested HF (high-frequency) sources, radio networks or the like can influence the operation of the device and may not be operated in combination with the system. Vario is a medical device that requires special safety precautions and must be installed and placed in operation in accordance with the attached EMC information. Portable and mobile RF communication devices (mobile telephones) can affect Vario.

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Electromagnetic Emissions		
Vario is intended for use in the electromagnetic environment specified below. The customer or the user of Vario should assure that they are used in such environment.		
Emission Tests	Compliance	Electromagnetic environment – guidance
RF emissions CISPR 11	Group 1	Vario uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference with nearby electronic equipment.
RF emissions CISPR 11	Class B	Vario is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Not applicable	

WARNINGS

Vario should not be used adjacent to or stacked with other equipment. If adjacent or stacked use is necessary, Vario should be observed to verify normal operation in the configuration in which it will be used.

Electromagnetic Immunity

Vario is intended for use in the electromagnetic environment specified below. The customer or the user of Vario should assure that they are used in such environment.

Immunity Tests	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Electrostatic Discharge (ESD) IEC 61000-4-2	± 6 kV contact ± 8 kV air	± 6 kV contact ± 8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrical fast transient/burst IEC 61000-4-4	± 2 kV for power supply lines ± 1 kV for input / output lines	± 2 kV for power supply lines not applicable	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	± 1 kV line(s) to line(s) not applicable	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5 % U_T (>95 % dip in U_T) for 0.5 cycle 40 % U_T (60 % dip in U_T) for 5 cycles 70 % U_T (30 % dip in U_T) for 25 cycles <5 % U_T (>95 % dip in U_T) for 5 sec	<5 % U_T (>95 % dip in U_T) for 0.5 cycle 40 % U_T (60 % dip in U_T) for 5 cycles 70 % U_T (30 % dip in U_T) for 25 cycles <5 % U_T (>95 % dip in U_T) for 5 sec	Mains power quality should be that of a typical commercial or hospital environment. If the user of Vario requires continued operation during power mains interruptions, it is recommended that Vario be powered from an uninterruptible power supply or a battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels of a typical commercial or hospital environment.

Note: U_T is the a.c. mains voltage prior to application of test level.

Electromagnetic Immunity

Vario is intended for use in the electromagnetic environment specified below. The customer or the user of Vario should assure that Vario is used in such environment.

Immunity Tests	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz	3 Vrms	Portable and mobile RF communications equipment should be used no closer to any part of Vario, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance $d = 1.2 \sqrt{P}$ $d = 1.2 \sqrt{P}$ 80 MHz to 800 MHz $d = 2.3 \sqrt{P}$ 800 MHz to 2.5 GHz where <i>P</i> is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and <i>d</i> is the recommended separation distance in meters (m) Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, should be less than the compliance level in each frequency range. ^b Interference may occur in the vicinity of equipment marked with the following symbol: 
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	3 V/m	

NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

^a Field strengths from fixed RF transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which Vario are used exceeds the applicable RF compliance level above, Vario should be observed to verify normal operation. If abnormal operation is observed, additional measures may be necessary, such as reorienting or relocating Vario.

^b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

Recommended separation distance between portable and mobile RF communications equipment and Vario

Vario is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of Vario can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and Vario as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter	Separation distance according to frequency of transmitter		
	150 kHz to 80 MHz $d = 1.2 \sqrt{P}$	80 MHz to 800 MHz $d = 1.2 \sqrt{P}$	800 MHz to 2.5 GHz $d = 2.3 \sqrt{P}$
W			
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

Notes

Notes

Medela AG
Lättichstrasse 4b
6341 Baar, Switzerland
www.medela.com

International Sales

Medela AG
Lättichstrasse 4b
6341 Baar
Switzerland
Phone +41 41 562 51 51
Fax +41 41 562 51 00
ism@medela.ch
www.medela.com

Australia

Medela Australia Pty Ltd,
Medical Technology
3 Arco Lane
Heatherton, Vic 3202
Australia
Phone +61 3 9552 8600
Fax +61 3 9552 8699
contact@medela.com.au
www.medela.com.au

Canada

Medela Canada Inc.
4160 Sladeview Crescent Unit # 8
Mississauga, Ontario
L5L 0A1
Canada
Phone +1 905 608 7272
Fax +1 905 608 8720
info@medela.ca
www.medela.ca

China

Medela (Beijing)
Medical Technology Co., Ltd.
Unit M, 4F
No 1068 Wuzhong Road,
Minhang District
Shanghai 201103
P.R. of China
Phone +86 21 62368102
Fax +86 21 62368103
info@medela.cn
www.medela.cn

France

Medela France Sarl
14, rue de la Butte Cordière
91154 Etampes cedex
France
Phone +33 1 69 16 10 30
Fax +33 1 69 16 10 32
info@medela.fr
www.medela.fr

Germany

Medela Medizintechnik
GmbH & Co. Handels KG
Postfach 1148
85378 Eching
Germany
Phone +49 89 31 97 59-0
Fax +49 89 31 97 59 99
info@medela.de
www.medela.de

India

Medela India private limited
c/o Vatika Business Park
First floor, tower 2,
Sohna Road, Sec-49,
Gurgaon 122 002
India
Phone +91 124 4416999
Fax +91 124 4416990
info@medela.in
www.medela.in

Italy

Medela Italia Srl
Via Turrini, 13/15 – Loc. Bargellino
40012 Calderara di Reno (BO)
Italy
Phone +39 051 72 76 88
Fax +39 051 72 76 89
info@medela.it
www.medela.it

Japan

Medela K.K.
KDX Kiba Building 3F
5-12-8, Kiba, Koto-ku,
Tokyo 135-0042
Japan
Phone +81 3 3820-5187
Fax +81 3 3820-5181
info@medela.jp
www.medela.jp

Netherlands & Belgium

Medela Benelux BV
Uilenwaard 31
5236 WB's-Hertogenbosch
Netherlands
Phone +31 73 690 40 40
Fax +31 73 690 40 44
info@medela.nl
info@medela.be
www.medela.nl
www.medela.be

Poland

Medela Polska Sp. z o.o.
Ul. Lewinowska 8
03-684 Warszawa
Poland
Phone +48 22 865 12 50
Fax +48 22 864 38 65
info@medela.pl
www.medela.pl

Russia

OOO Medela
Vavilova str. 97, entrance 3
117335 Moscow
Russia
Phone +7 495 980 6194
Fax +7 495 980 6194
info@medela-russia.ru
www.medela-russia.ru

Spain & Portugal

Productos Medicinales Medela, S.L.
c/ Manuel Fernández Márquez, 49
08918 Badalona (Barcelona)
Spain
Phone +34 93 320 59 69
Fax +34 93 320 55 31
info@medela.es
info@medela.pt
www.medela.es
www.medela.pt

Sweden

Medela Medical AB
Box 7266
187 14 Täby
Sweden
Phone +46 8 588 03 200
Fax +46 8 588 03 299
info@medela.se
www.medela.se

Switzerland

Medela AG
Lättichstrasse 4b
6341 Baar
Switzerland
Phone +41 41 562 51 51
Fax +41 41 562 51 00
contact@medela.ch
www.medela.ch

United Kingdom

Medela UK Ltd.
Huntsman Drive
Northbank Industrial Park
Irlam, Manchester M44 5EG
UK
Phone +44 870 950 5994
Fax +44 870 389 2233
info@medela.co.uk
www.medela.co.uk

USA

Medela Inc.
1101 Corporate Drive
McHenry, IL 60050
USA
Phone +1 877 735 1626
Fax +1 815 307 8942
suction@medela.com
www.medelasuction.com