EN

DV6WM DeVilbiss Wireless Modem Instruction Guide

Made in USA of US and Imported Parts.

ES

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INTENDED USE

The DeVilbiss Healthcare DV6WM Wireless Modem is intended to be used as a data collection tool for DeVilbiss IntelliPAP®/SleepCube® Series with SmartLink® II Bluetooth® module, IntelliPAP 2, and DeVilbiss BLUE™ Series CPAPs for patients in a home or healthcare environment. It is not intended to be used as a diagnostic tool.




CONTRAINDICATIONS

 MR Unsafe

Do not bring the device or accessories into a Magnetic Resonance (MR) environment as it may cause unacceptable risk to the patient or damage to the device or MR medical devices. The device and accessories have not been evaluated for safety in an MR environment.

Do not use the device or accessories in an environment with electromagnetic equipment such as CT scanners, Diathermy, RFID and electromagnetic security systems (metal detectors) as it may cause unacceptable risk to the patient or damage to the device. Some electromagnetic sources may not be apparent, if you notice any unexplained changes in the performance of this device, if it is making unusual or harsh sounds, disconnect the power cord and discontinue use. Contact your home care provider.

SYMBOLS

	Test Call Button		Manufacturer
	Cell Phone Status LED		Date of Manufacture
	Bluetooth Status LED		DC voltage symbol
	Battery Status LED	REF	Model number
	Battery Charge Connector	SN	Serial Number
	Refer to Operating Instructions	EC REP	European Representative
	Radio transmitter		MR Unsafe – Unsafe for Magnetic Resonance Environment
	IP21	Ingress Protection - Protected against finger access to hazardous parts; protected against vertically falling water drops	
Rx Only	CAUTION - USA Federal law restricts this device to sale by or on the order of physician		
	This device contains electrical and/or electronic equipment that must be recycled per EU Directive 2012/19/EU- Waste Electrical and Electronic Equipment (WEEE)		

IMPORTANT SAFEGUARDS

READ ALL INSTRUCTIONS BEFORE USING THIS DEVICE.

SAVE THESE INSTRUCTIONS.

(per IEC 60601-1-2:2007 Section 5.2.2.1 b)

WARNING: The use of ACCESSORIES, transducers and cables other than those specified, with the exception of transducers and cables sold by DeVilbiss Healthcare as replacement parts for internal components, may result in increased EMISSIONS or decreased IMMUNITY of the ME EQUIPMENT or ME SYSTEM.

(per IEC 60601-1-2:2007 Section 5.2.2.1 d)

WARNING: The DeVilbiss DV6WM Wireless Modem should not be used adjacent to or stacked with other equipment and that if adjacent or stacked use is necessary, the DeVilbiss DV6WM Wireless Modem should be observed to verify normal operation in the configuration in which it will be used.

(per IEC 60601-1-2:2007 Section 5.2.2.5)

WARNING: The DeVilbiss DV6WM Wireless Modem may be interfered with by other equipment, even if that other equipment complies with CISPR EMISSION requirements.

RF receiver frequency band BT	2.402 GHz ~ 2.480 GHz
RF receiver bandwidth BT	2.402 GHz ~ 2.480 GHz
RF receiver frequency band GSM	824.2 MHz ~ 1907.6 MHz
RF receiver bandwidth GSM	824.2 MHz ~ 1907.6 MHz

(per IEC 60601-1-2:2007 Section 5.2.2.5)

RF transmit frequency band BT	2.402 GHz ~ 2.480 GHz
Modulation BT	GFSK, 8DPSK
Effective Radiated Power BT	< -10 dBm (100 mW)
RF transmit frequency band GSM	824.2 MHz ~ 1907.6 MHz
Modulation GSM	8-PSK, WCDMA
Effective Radiated Power GSM	< +33 dBm(2,000 mW)

CAUTIONS / WARNINGS

- Refer to International Standard IEC 60601-1 Ed 3.0 Amendment 1 for safety requirements applicable to Medical Electrical Systems.
- Electric shock hazard-Do not attempt to open or remove the Wireless Modem or charger cabinet; there are no user-serviceable components inside. If service is required, contact your equipment provider for instructions on obtaining service. Opening or attempting to service your device will void the warranty.
- Do not perform service or maintenance while the Wireless Modem is in use.
- Wireless Modem contains Lithium Polymer (LiPO) battery. Do not dispose of in household waste.
- The replacement of lithium batteries by inadequately trained personnel could result in a HAZARD. There are no user serviceable parts inside the DV6WM Wireless Modem, the modem should only be serviced by a qualified DeVilbiss provider.
- If the Wireless Modem begins to perform in an unusual manner or a way not described in this guide, do not attach the Wireless Modem to a charger or any other equipment. Contact your equipment provider and report the unusual performance.
- If Wireless Modem is exposed to extreme temperatures or humidity during storage, allow the device to sit at room temperature for 2 hours before operating or charging.
- Do not connect the DV6WM Wireless Modem to equipment not described in these instructions. Connecting the modem to other equipment or accessories may result in malfunction or damage to the DV6WM Wireless Modem.
- Do not modify the DV6WM Wireless Modem or accessories this may result in malfunction or damage to the DV6WM Wireless Modem.

CAUTIONS / WARNINGS

- Use only chargers and charging cables provided by DeVilbiss for use with the DV6WM Wireless Modem. Use of other chargers or charging cables may result in malfunction or damage to the DV6WM Wireless Modem.

For storage and handling – The following conditions may cause a fire hazard by causing a failure of the internal rechargeable battery:

- Never attempt to charge a Wireless Modem with physical damage.
- If Wireless Modem temperature exceeds 50°C (122°F), terminate charging.
- Do not heat or throw Wireless Modem into a fire.
- Do not use or store Wireless Modem close to fire or inside the car where temperature may be over 60°C.
- Do not immerse the Wireless Modem in water.

PACKAGE CONTENTS

DV6WM Wireless Modem


USB Wall Charger

Instruction Guide

IMPORTANT PARTS

Wireless Modem – Control Panel

Normal Operation

The battery LED flashes when battery is low, otherwise all LEDs remain OFF except during call sequence. During the call, the battery LED illuminates solid and either Bluetooth or Cell LEDs are also active. You may not see this happen.

Call frequency has been preset by your provider.


NOTE: When battery LED flashes, refer to *Charging the Modem Battery*.

NOTE: To check operation at any time, press the test call button and confirm LEDs illuminate.


SET UP & OPERATION - DV6x Series CPAP

Initial Set Up of Your Wireless Modem

NOTE: Fully charge battery before beginning set up. Refer to **Charging the Modem Battery**.



1. Check DV6X CPAP to ensure Bluetooth is ON
 - a. Navigate to ADD-ONS
 - b. Select **Bluetooth**
 - c. Press to toggle from OFF to ON if necessary






2. Press **Test Call** button on Wireless Modem

- a. All LEDs will illuminate for approximately 1 second, then turn OFF
- b. Battery LED shows battery status
- c. Bluetooth LED will begin flashing slowly while searching for CPAP

3. Pair with DV6x Series CPAP
 - a. Modem bluetooth LED will begin flashing fast when ready to pair with CPAP (approximately 30 seconds for fast flash to begin)
 - b. The wireless modem serial number will appear on the list when the CPAP finds it

NOTE: Wait until Modem appears in list, this may take 30 seconds.

- c. At Click Item to Pair screen, select **Wireless Modem** serial number and press, the menu will show "Pairing..." then "Success"
- d. The Bluetooth LED on the Wireless Modem will stop flashing and stay ON while reading data from the CPAP

4. The Cell Phone LED on the modem will flash and then turn solid for a few minutes


5. Find a location where Cell Phone LED is ON (not flashing) and place the modem where it will not be disturbed

NOTE: Bluetooth and Cell Phone LEDs will normally be OFF; they will be active only when modem is talking to CPAP or to the server; you may not see this happen. Leave modem alone.


6. After call sequence is finished, all LEDs will turn OFF

SET UP & OPERATION - DV6x Series CPAP

Operating your Wireless Modem with Your DV6x Series CPAP


The DeVilbiss Wireless Modem operates without interaction as long as the battery level is good.


To see modem info, select **Wireless Modem** in the ADD ONS menu on the DV6x CPAP and press.

Info shown on DV6x CPAP display for Wireless Modem


- Days since last upload to server
- Wireless Modem ID
- Bluetooth Signal
- Cell Phone Signal
- Battery Level


NOTE: The DV6x Series CPAP will report low modem battery. When you see this message, refer to **Charging the Modem Battery**. The low battery message appears when blower is turned ON.

CHARGING THE MODEM BATTERY

When the battery level gets low, the Battery LED on the modem will flash. Refer to **Charging the Modem Battery**.

The Battery LED will flash 3 times in 1 second every 3 seconds when the battery level gets low. If the Battery LED is flashing, plug the provided USB Wall Charger into the modem and the wall outlet.

Battery LED on the modem will:

- Flash slowly when battery is charging
- Stay on steady when battery is fully charged and connected to the power source

NOTE: *If desired, the modem can be plugged in at all times.*

The DV6WM Wireless modem operates normally while the battery is charging. The modem may be left plugged into the charger for an extended period of time.

Typical operation time of the DV6WM Wireless Modem on a full battery charge is 90 days.

CLEANING INSTRUCTIONS

1. Clean modem once per month. Unplug Wireless Modem from charger and wipe outside surface with a damp cloth.
2. Store Wireless Modem in a dry place at room temperature when not in use.

Multi-Patient Use

If the device is to be used by multiple persons in a healthcare environment or device rental program, the outside of the device must be cleaned and disinfected between patients.

Cleaning for Multi-Patient Use

WARNING

If you are using the device on multiple patients, complete the following steps to clean the outside of the device only before each new user.

1. Unplug the device before cleaning.
2. Use a cloth with a solution of mild detergent (e.g. Dawn dishwashing liquid) and water to remove debris if necessary.
3. Use a cloth with a 1:10 dilution of 5.25%–6.15% sodium hypochlorite (i.e., household bleach) to clean the exterior of the device.
4. Ensure that the device is completely dry before plugging in the power cord.

MAINTENANCE

WARNING

Electric shock hazard-Do not attempt to open or remove the Wireless Modem or charger cabinet; there are no user-serviceable components inside. If service is required, contact your equipment provider for instructions on obtaining service. Opening or attempting to service your device will void the warranty.

The Wireless Modem and accessories require no maintenance through the expected service life.

EXPECTED SERVICE LIFE

DV6WM Wireless Modem	5 years
AC to DC USB wall charger	5 years

Battery Replacement – The battery in the DV6WM Wireless Modem is not user replaceable. If your battery needs to be replaced, contact your equipment provider for service.

DISPOSAL

CAUTION: *Wireless Modem contains Lithium Polymer (LiPO) battery. Do not dispose of in household waste.*

Wireless Modem and charger contain electronic waste, contact local authorities to determine proper method of disposal of electronic parts, plastic parts and LiPO battery.

NETWORK/DATA COUPLING

- Connection of the DV6WM Wireless Modem to a network/data coupling that includes other equipment could result in previously unidentified risks to patients and operators, and the responsible organization should identify, analyze, and control such risks;
- Subsequent changes to network/data coupling may introduce new risks, requiring new analysis;

Changes to network/data coupling include:

- Changes in network/data coupling configuration
- Connection of additional items to network/data coupling
- Disconnecting items from network/data coupling
- Update of equipment connected to network/data coupling
- Upgrade of equipment connected to network/data coupling

SPECIFICATIONS

Size	4.4"H x 3"W x 1"D (11.2cm x 7.6cm x 2.5cm)
Weight	0.35 lbs. (158.8 g)
Electrical Requirements AC	100 – 240 VAC, 50/60 Hz, 120 mA (Specified external AC to DC USB charger input)
Electrical Requirements DC	4.74 to 5.25 VDC, 500 mA (Specified external AC to DC USB charger output)
Maximum Power Consumption	3 watts from AC to DC USB charger
Operating Temperature Range	41°F to 104°F (5°C to 40°C)
Operating Humidity Range	15% to 95% RH non-condensing
Operating Atmospheric Conditions	700 hPa to 1060 hPa (~9800 ft to ~ 1400ft below sea level)
Storage & Transportation Temperature Range	13°F to +158°F (-25°C to +70°C)
Storage & Transportation Humidity Range	0% to 95% RH non-condensing

Technical Specs (per IEC 60601-1-2:2007 Section 5.2.2.1 a)

Cables compatible with the DeVilbiss DV6WM Wireless Modem.....	USB to Micro B USB cable, max cable length = 6 ft (1.8m)
USB Wall Charger	Input 100 - 240 VAC, 47 – 63 Hz Output 5 VDC, 3 Watt max Output connector USB Micro B Max cable length = 6 ft (1.8 m)

Wireless Specifications: Bluetooth

This medical device contains two radio transmitters, Bluetooth and 3G cellular. The Bluetooth radio integrated into this device is active when the Bluetooth Status LED (⌘) on the control panel is illuminated. Bluetooth is used to connect wirelessly to your medical device. The 3G cellular radio integrated into this device is active when the Cell Phone Status LED (☎) on the control panel is illuminated. The 3G cellular radio is used to send data read from your medical device to a remote data server. If you notice any unexplained changes in the performance of the wireless function or your device, turn the Bluetooth radio off to see if that is the cause. You may experience issues with wireless technology such as delays in displayed or stored data, artifact or gaps in the data due to information lost during wireless transfer.

Radio Technology.....	Bluetooth 2.1 +EDR and Bluetooth 4.0
Bluetooth Power Class	1.5
Network Topology.....	Point to Point
Bluetooth Profile Supported	SPP
Effective Range	50 meters (line of sight)
Effective Radiated Power	10 dBm (100 mW)
Radio frequency band (Tx and Rx).....	2.402 GHz ~ 2.480 GHz
Minimum Separation Distance (to other RF transmitters).....	1 cm (0.4 in.)

Wireless Specifications 3G Cellular:

Radio Technology	3G Cellular
Network Topology	Point to Point
Effective Radiated Power	max 33 dBm (2 W)
Radio frequency band (Tx and Rx)	824.2 MHz ~ 1907.6 MHz
Minimum Separation Distance (to other RF transmitters)	1 cm (0.4 in.)
Quality Of Service Required	1 Bar (> -124 dBm)

Security Requirements:

Authentication	Enforced on all data channels (outgoing and incoming)
Encryption	Base 64 encoding

TROUBLESHOOTING

CAUTION: If the Wireless Modem begins to perform in an unusual manner or a way not described in this guide, do not attach the Wireless Modem to a charger or any other equipment. Contact your equipment provider and report the unusual performance.

CAUTION: If Wireless Modem is exposed to extreme temperatures or humidity during storage, allow the device to sit at room temperature for 2 hours before operating or charging.

Issue	Possible Cause	Remedy
Nothing happens after pressing the test call button.	1. Wireless Modem battery is discharged.	1. Plug Wireless Modem into charger and allow battery to charge for at least 2 hours.
	2. Test call button is not working properly.	2. & 3. Contact your equipment provider for service.
	3. Malfunction of the modem.	
Unable to pair Wireless Modem with a CPAP.	1. CPAP does not have Bluetooth.	1. Wireless Modem only works with DeVilbiss CPAPs that have Bluetooth feature.
	2. CPAP Bluetooth is not turned on.	2. Turn Bluetooth on in CPAP menu, see instructions.
	3. Wireless Modem is too far away from CPAP.	3. Move modem closer to CPAP, must be within 30m of each other.
Battery LED flashes fast three times in a row.	1. Battery charge is low.	1. Allow Wireless Modem to sit at room temperature for 2 hours, <u>attach charger to charge battery.</u>
Wireless Modem does not charge, Battery LED does not light when connected to charger.	1. Charger may not be working	1. Try using another charger, one that is approved for use with the Wireless Modem.
	2. Malfunction of Wireless Modem	2. Contact your equipment provider for service.
LEDs flashing in a pattern not explained in this guide.	1. Malfunction of Wireless Modem	1. Contact your equipment provider for service.

Contact your equipment provider if you need assistance setting up, using or maintaining your Wireless Modem.

Report unexpected operation or events to your equipment provider or DeVilbiss Customer Service 800-338-1988. See back cover of guide for other DeVilbiss office locations.

EMC INFORMATION

(per IEC 60601-1-2:2007 Section 5.2.1.1)

MEDICAL ELECTRICAL EQUIPMENT needs special precautions regarding EMC and needs to be installed and put into service according to the EMC information provided in this instruction guide.

Portable and mobile RF communications equipment can affect MEDICAL ELECTRICAL EQUIPMENT

(per IEC 60601-1-2:2007 Section 5.2.2.1 c)

Guidance and Manufacturer's Declaration – Electromagnetic Emissions		
The DeVilbiss DV6WM Wireless Modem is intended for use in the electromagnetic environment specified below. The customer or the user of the DeVilbiss DV6WM Wireless Modem should assure that it is used in such an environment.		
Emissions Test	Compliance	Electromagnetic Environment – Guidance
RF Emissions CISPR 11	Group 2	This device uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF Emissions CISPR 11	Class B	This device is suitable for use in all establishments including domestic, and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonics IEC 61000-3-2	Class A	
Voltage Fluctuations/ Flicker Emissions IEC 61000-3-3	Complies	

(per IEC 60601-1-2:2007 Section 5.2.2.1 f)

Immunity Test	IEC 60601 test level	Compliance level	Electromagnetic Environment – Guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±8kV contact ±15kV air	±8kV contact ±15kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrical fast transient/burst IEC 61000-4-4	±2 kV for power supply lines ±1 kV for input/output lines	±2 kV for power supply lines ±1 kV for input/output lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	±1 kV line(s) to line(s) ±2 kV line(s) to earth	±1 kV line(s) to line(s) ±2 kV line(s) to earth	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	>95 % dip in UT for 0.5 cycle 60 % dip in UT for 5 cycles 30 % dip in UT for 25 cycles >95 % dip in UT for 5 seconds	>95 % dip in UT for 0.5 cycle 60 % dip in UT for 5 cycles 30 % dip in UT for 25 cycles >95 % dip in UT for 5 seconds	Mains power quality should be that of a typical commercial or hospital environment. If the user of the DeVilbiss DV6WM Wireless Modem requires continued operation during power mains interruptions, it is recommended that DeVilbiss DV6WM Wireless Modem be powered from an uninterruptible power supply or a battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	30 A/m	30 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

NOTE: *UT is the a.c. mains voltage prior to application of the test level.*

EMC INFORMATION

(per IEC 60601-1-2:2007 Section 5.2.2.2)

Guidance and Manufacturer's Declaration – Electromagnetic Immunity			
The DeVilbiss DV6WM Wireless Modem is intended for use in the electromagnetic environment specified below. The customer or the user of the DeVilbiss DV6WM Wireless Modem should assure that it is used in such an environment.			
Immunity Test	IEC 60601 test level	Compliance level	Electromagnetic Environment – Guidance
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz	V1 = 3 Vrms 6 Vrms on ISM & Amateur Bands	Portable and mobile RF communications equipment should be separated from the device by no less than the recommended separation distances calculated/listed below: $D=(0.4)\sqrt{P}$
Radiated RF IEC 61000-4-3	10 V/m 80 MHz to 2.5 GHz	E1 = 10 V/m	$D=(0.4)\sqrt{P}$ 80 to 800 MHz $D=(0.7)\sqrt{P}$ 800 MHz to 2.5 GHz Where P is the maximum power rating in watts and D is the recommended separation distance in meters. Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey ^a , should be less than the compliance level in each frequency range ^b . Interference may occur in the vicinity of equipment marked with the following symbol: 
<p>NOTE 1: At 80 MHz and 800 MHz, the higher frequency range applies.</p> <p>NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.</p> <p>a. Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the DeVilbiss DV6WM Wireless Modem is used exceeds the applicable RF compliance level above, the DeVilbiss DV6WM Wireless Modem should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the DeVilbiss DV6WM Wireless Modem.</p> <p>b. Over the frequency range 150 kHz to 80 MHz, field strengths should be less than [V1] V/m.</p>			

(per IEC 60601-1-2:2007 Section 5.2.2.2)

Recommended Separation Distances Between Portable and Mobile RF Communications Equipment and this device. This device and system are NOT Life-Supporting			
The DeVilbiss DV6WM Wireless Modem is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the DeVilbiss DV6WM Wireless Modem can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the DeVilbiss DV6WM Wireless Modem as recommended below, according to the maximum output power of the communications equipment.			
Rated Maximum Output Power of Transmitter (W)	Separation distance according to frequency of transmitter (m)		
	150 kHz to 80 MHz outside ISM bands $D=(0.4)\sqrt{P}$	80 MHz to 800 MHz $D=(0.4)\sqrt{P}$	800 MHz to 2.5 GHz $D=(0.7)\sqrt{P}$
0.01	0.04	0.04	0.07
0.1	0.11	0.11	0.22
1	0.35	0.35	0.70
10	1.1	1.1	2.2
100	3.5	3.5	7.0
For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be determined using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.			
<p>NOTE 1: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.</p> <p>NOTE 2: The ISM (industrial, scientific and medical) bands between 150 kHz and 80 MHz are 6.765 MHz to 6.795 MHz; 13.553 MHz to 13.567 MHz; 26.957 MHz to 27.283 MHz; and 40.66 MHz to 40.70 MHz.</p> <p>NOTE 3: An additional factor of 10/3 has been incorporated into the formulae used in calculating the recommended separation distance for transmitters in the ISM frequency bands between 150 kHz and 80 MHz and in the frequency range 80 MHz to 2.5 GHz to decrease the likelihood that mobile/portable communications equipment could cause interference if it is inadvertently brought into patient areas.</p> <p>NOTE 4: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.</p>			

EMC INFORMATION

FCC and Industry Canada Statement

GSM Module Contains FCC ID: QIPEHS6; Contains IC: 7830A-EHS6

Bluetooth Module Contains FCC ID: TFB-BT1; Contains IC: 5969A-BT1

This device complies with part 15 of the FCC Rules and with Industry Canada license-exempt RSS standard(s). Operation is subject to the following two conditions: (1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

The output power of this device is less than 20mW. The SAR test is not required. When using it, ensure that the antenna of the device is at least 20cm away from all persons.

This Class B digital apparatus complies with Canadian ICES-003.

NOTE: *The manufacturer is not responsible for any radio or TV interference caused by unauthorized modifications or changes to this equipment. Such modifications or changes could void the user's authority to operate the equipment.*

NOTE: *This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses, and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment causes harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:*

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

CE Declaration

Hereby, DeVilbiss Healthcare LLC, declares that this DV6WM Wireless Modem is in compliance with the essential requirements and other relevant provisions of Directive 1999/5/EC. The declaration of conformity may be consulted at www.devilbisshealthcare.com/files/DV6WM_CE_Declaration.pdf.

The Bluetooth® word, mark, and logos are registered trademarks owned by Bluetooth SIG, Inc. and any use of such marks by DeVilbiss is under license. Other trademarks and trade names are those of their respective owners.

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DeVilbiss®, IntelliPAP®, SleepCube®, and SmartLink® are registered trademarks of DeVilbiss Healthcare. Bluetooth® is a registered trademark of Bluetooth SIG, Inc. DeVilbiss BLUE™ is a trademark of DeVilbiss Healthcare.
© 2016 DeVilbiss Healthcare LLC. 03.16 All Rights Reserved.
DeVilbiss Healthcare Ltd
DeVilbiss Healthcare Pty. Limited
DeVilbiss Healthcare SAS
DeVilbiss Healthcare GmbH